

NEW!

BufferPro™

8.4% Sodium Bicarbonate
Buffer Solution



Faster, Easier, Better! One-Step, Simple Buffering!



BufferPro™ offers an elegant, one-step, affordable solution for buffering your local anesthetic injections. BufferPro dispenses 0.1mL of sodium bicarbonate into a dental anesthetic cartridge to raise the pH of the anesthetic solution close to physiologic levels. BufferPro is a single-use, sterile packaged, self-contained capsule with no measuring, mixing, or additional accessories.

#01B0000: BufferPro™ Solution, 25 Capsule/Box



Predictably Fast Buffering
One-step, simple buffering



Single-Use, Sterile Packaged
*Avoids cross-contamination,
preserves shelf-life*



Better Experience
*No accessories, mixing,
or measuring*



Economically Priced
Low cost per injection



septodontusa.com
from the manufacturer of Septocaine®